

[Medical and Surgical Supplies](#)

[Medical Disposables](#)

[Medical Equipment and Furniture](#)

[Medical Equipment Hire : Home Based Care](#)

Vat No.4330218381 CK 2004/114486/23

Physical: 138 Keal Road Sydenham Durban KZN SA 4091

Postal P O Box 37908 Overport Durban KZN SA 4067

Tel : 031 2085831/031 207 8354 Fax : 031 208 5844/0866 01 9791

E-mail : [info@alphahealthcare.co.za](mailto:info@alphahealthcare.co.za) Website : [www.alphahealthcare.co.za](http://www.alphahealthcare.co.za)

## DIAGNOSTICS RAPID TESTS

**Rapid In Vitro Diagnostics (IVD) test kit supplier.**

**We sell high quality hCG-Pregnancy, HIV, Syphilis, and Multiple Drug test kits, which are manufactured in South Africa under strict ISO 13485:2003 conditions.**

### **HCG - Pregnancy Test Kit**

Item Number: AHCE**540001**

**One Step / 3 Minute / Early Detection / Highly Accurate**  
Rapid Test kit used for qualitative determination of human chorionic gonadotropin (HCG) in urine for early detection of pregnancy.



### **HIV - HIV 1&2 Test Kit**

- For the rapid detection of HIV 1 & 2 antibodies in serum & whole blood
- Highly accurate results within 15 minutes

Item Number: AHCE**540002**



### **MDT - Multiple Drug Test Kit**

Item Number: AHCE**540003**

**Urine Based / One Step / Highly Accurate / Results in 5 Minutes**

Can detect drug metabolites in urine up to 72 hours after drug usage, at the following cut-off levels:

- AMP (d-Amphetamine) 1000 ng/ml → *Tik (crystal meth, glass, speed, ice)*
- COC (Benzoyllecgonine) 300 ng/ml → *Cocaine (crack, coke, rocks, snow)*
- THC (Marijuana) 50 ng/ml → *Marijuana (dagga, weed, pot, dope, hashish)*





Vat No.4330218381 CK 2004/114486/23

Physical: 138 Keal Road Sydenham Durban KZN SA 4091

Postal P O Box 37908 Overport Durban KZN SA 4067

Tel : 031 2085831/031 207 8354 Fax : 031 208 5844/0866 01 9791

[Medical and Surgical Supplies](#)

[Medical Disposables](#)

[Medical Equipment and Furniture](#)

[Medical Equipment Hire : Home Based Care](#)

E-mail : [info@alphahealthcare.co.za](mailto:info@alphahealthcare.co.za) Website : [www.alphahealthcare.co.za](http://www.alphahealthcare.co.za)

**Syphilis -Syphilis Test Kit**

Item Number: AHCE**540004**



**You may request a Quote and Specifications**