

Medical and Surgical Supplies

Medical Disposables

Medical Equipment and Furniture

Medical Equipment Hire : Home Based Care

E-mail : info@alphahealthcare.co.za Website : www.alphahealthcare.co.za

Vat No.4330218381 CK 2004/114486/23

Physical: 138 Keal Road Sydenham Durban KZN SA 4091

Postal P O Box 37908 Overport Durban KZN SA 4067

Tel : 031 2085831/031 207 8354 Fax : 031 208 5844/0866 01 9791



Littmann Classic II S.E. Stethoscope Model 2201

Product Description: Colour : Black

The 'industry standard' stethoscope. It features 3M's patented tuneable diaphragm which enables you to listen to high and low frequency sounds by simply alternating pressure on the chest piece.

- excellent acoustic seal and comfortable fit with patented Snap Tight Soft-Sealing Ear tips
- non-chill rims for patient comfort
- comfortably angled, anatomically correct headset
- high performance and extremely long-lasting single-lumen tubing
- stainless steel chest pieces
- 71cm total length
- available in a variety of colours
- 3 year warranty

Promotional Price – R1190,00 incl Vat



Add Name Engraving (optional) for R55.00