



[Medical and Surgical Supplies](#)

[Medical Disposables](#)

[Medical Equipment and Furniture](#)

[Medical Equipment Hire : Home Based Care](#)

**E-mail :** [info@alphahealthcare.co.za](mailto:info@alphahealthcare.co.za) **Website :** [www.alphahealthcare.co.za](http://www.alphahealthcare.co.za)

**Vat No.4330218381 CK 2004/114486/23**

**Physical:** 138 Keal Road Sydenham Durban KZN SA 4091

**Postal** P O Box 37908 Overport Durban KZN SA 4067

**Tel :** 031 2085831/031 207 8354 **Fax :** 031 208 5844/0866 01 9791

## **Beds**

### **Electric / Gas Lift / Winding**

### **Home Based Care**



#### **Long-Term Care Bed**

**Long-Term Care (LTC)** bed offers a comfortable working height for caregivers, while helping to facilitate easy resident ingress and egress.

The Resident® LTC bed provides a unique combination of function, aesthetics, durability, and value to enhance the well being of both resident and facility. Separate articulation of the head and knee areas helps to personalize comfort.

**Electric High/Low Controls** – These easy-to-use controls facilitate resident egress and provide a comfortable working height for staff.

**Half-Length Tuck Away Side Rails** – Side rails fold out of the way of knees and shins.

**Easy-to-use and see Operating Controls** – These controls contribute to resident comfort and independence.

**Automatic Contour** – This feature helps protect fragile skin and facilitates resident management.

**Foot end Control Lockout Panel** – May be used to deny resident access to operating functions.

**Floor Brake** – Helps to prevent unintended bed movement that can pose risk to patient safety.

Medical and Surgical Supplies

Medical Disposables

Medical Equipment and Furniture

Medical Equipment Hire : Home Based Care

Vat No.4330218381 CK 2004/114486/23

Physical: 138 Keal Road Sydenham Durban KZN SA 4091

Postal P O Box 37908 Overport Durban KZN SA 4067

Tel : 031 2085831/031 207 8354 Fax : 031 208 5844/0866 01 9791

E-mail : [info@alphahealthcare.co.za](mailto:info@alphahealthcare.co.za) Website : [www.alphahealthcare.co.za](http://www.alphahealthcare.co.za)

**M1**

**HI LO BED - GASLIFT BACKREST AND TRENDELENBURG**



Epoxy coated construction with removable head and foot ends. Backrest adjustable by means of gas lift lever allocated on the backrest of the bed. Gas spring assisted forward and reverse trendelenburg Collapsible safety sides.

**M1A**

**HI-LO BED WINDING BACKREST AND TRENDELENBURG**



Epoxy coated construction with removable head and foot ends. Backrest adjustable by means of winding handle allocated by the footed of the bed. Gas spring assisted forward and reverse trendelenburg Collapsible safety sides.